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

Artificial skin Penetration test (Parafilm® M)


Insertion depths of NEOPATH™ (600 µm) in Parafilm® M (a) and in neonatal porcine skin (b) applied by a Texture Analyser (10 N).

	inserti
	Parafilm® M
Average	436.54
SD	11.58
RSD	2.65 %


6.22%

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

Microplastic Detection Test Based on Usage Frequency


Before insertion


After 100 insertions (10 N) on Parafilm®

Punctured parafilm® M

	Weight (mg)	
	Before insertion	After 100 insertions
Average	211.54	211.54
SD	0.02	0.01
RSD	0.01%	0.01%

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



Anjani, Qonita Kurnia, et al. "Parafilm® M and Strat-M® as skin simulants in in vitro permeation of dissolving microarray patches loaded with proteins." *International Journal of Pharmaceutics* 655 (2024): 124071.

Weimer, Patrícia, Rochele Cassanta Rossi, and Leticia Scherer Koester. "Dissolving microneedles developed in association with Nanosystems: A scoping review on the quality parameters of these emerging Systems for Drug or protein transdermal delivery." *Pharmaceutics* 13.10 (2021): 1601.

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
Larrañeta, E., Moore, J., Vicente-Pérez, E. M., González-Vázquez, P., Lutton, R., Woolfson, A. D., & Donnelly, R. F. A proposed model membrane and test method for microneedle.

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


NanoPass Technologies

Growth Score **86** ↓ CB Rank **61676** ↓ Heat Score **65** ↓

NanoPass is continuously engaged in joint clinical development with top pharmaceutical companies.

Founded 2000
 Private
 Seed
 Ness Ziona, HaMerkaz, Israel
 1-10
 www.nanopass.com/

[Clinical Trials](#)
[Medical](#)
[Pharmaceutical](#)

Financials

Company Funding

Funding Rounds

Investors

Total Funding Amount
\$10.5M

Number of Funding Rounds
2

Lead Investors
Awz Ventures, Ofer Hi-Tech, WFD Ventures, D Capital Partners

Announced Date	Transaction Name	Number of Investors	Money Raised	Lead Investors
Jul 15, 2024	Seed Round - NanoPass Technologies	1	\$4M	Awz Ventures
Feb 18, 2007	Seed Round - NanoPass Technologies	3	\$6.5M	—

[View Financial Details >](#)

Company Funding

Funding Rounds

Investors

Investor Name	Lead Investor	Funding Round	Partners
Awz Ventures	Yes	Seed Round - NanoPass Technologies	Yaron Ashkenazi
Ofer Hi-Tech	No	Seed Round - NanoPass Technologies	—
WFD Ventures	No	Seed Round - NanoPass Technologies	—
D Capital Partners	No	Seed Round - NanoPass Technologies	—

[View Financial Details >](#)

Predictions & Insights

Growth Insight ⓘ
updated Dec 27, 2024

Growing

NanoPass Technologies is very likely to be growing.

NanoPass Technologies secured a \$4 million investment from Awz Ventures, indicating growth and innovation in their intradermal delivery technology.

Top Contributing Factors

Recently Raised Funding >

[View Details >](#)

Growth Prediction ⓘ
updated Jan 10, 2025

95 Probable Future Growth

Top Contributing Factors

Runway Speculation

[View Details >](#)

IPO Prediction ⓘ
updated Feb 6, 2025

Doubtful

Top Contributing Factors

Growth Insight >

Acquisition Prediction ⓘ
updated Feb 6, 2025

Uncertain

Top Contributing Factors

Growth Insight >

Runway Speculation

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INFLUENZA

Van Damme P, et al. Safety and efficacy of a novel microneedle device for dose sparing intradermal influenza vaccination in healthy adults. *Vaccine* 2009; 27:454–9;

Della Cioppa G, et al. A dose–range study in older adults to compare the safety and immunogenicity profiles of MF59 adjuvanted and non–adjuvanted seasonal influenza vaccines following intradermal and intramuscular administration. *Hum Vaccines Immunother* 2014; 10:6, 1–10. Read full article.

Levin Y, Kochba E, Kenney R. Clinical evaluation of a novel microneedle device for intradermal delivery of an influenza vaccine: Are all delivery methods the same? *Vaccine* 2014; 34:4249–52;

Hung IF et al . Quantitative and qualitative analysis of antibody response after dose sparing intradermal 2009 H1N1 vaccination. *Vaccine* 30 (2012) 2707– 2708. Read full article.

Hung IF, et al. Dose sparing intradermal trivalent influenza (2010/2011) vaccination overcomes reduced immunogenicity of the 2009 H1N1 strain. *Vaccine* 2012; 30:6427–35;

Levin Y. et al. A phase 1, open–label, randomized study to compare the immunogenicity and safety of different administration routes and doses of virosomal influenza vaccine in elderly. *Vaccine* 34 (2016) 5262–5272. Read full article.

Carter et al. The adjuvant GLA–AF enhances human intradermal vaccine responses. *Sci. Adv.* 2018; 4. eaas9930. Read full article.

El Sahly HM, et al. Topical Imiquimod Does Not Provide an Adjuvant Effect When Administered with Inactivated Influenza A/H5N1 Vaccine in Healthy Young Adults, *The Journal of Infectious Diseases*, Volume 224, Issue 10, 15 November 2021, Pages 1712–1719.

POLIO

Stephanie B. Troy et al. Comparison of the Immunogenicity of Various Booster Doses of Inactivated Polio Vaccine Delivered Intradermally Versus Intramuscularly to HIV–Infected Adults. *J Infect Dis.* (2015) 211 (12): 1969–1976. Read full article..

Anand A et al. Early priming with inactivated poliovirus vaccine (IPV) and intradermal fractional dose IPV administered by a microneedle device: A randomized controlled trial. *Vaccine* 33 (2015) 6816–6822

ZOSTERC

han R Beals et al. Immune response and reactogenicity of intradermal administration versus subcutaneous administration of varicella–zoster virus vaccine: an exploratory, randomised, partly blinded trial. *Lancet Infect Dis.* 2016 Aug;16(8):915–22

HBV

Hung IF et al. A Double–blind, Randomized Phase 2 Controlled Trial of Intradermal Hepatitis B Vaccination With a Topical Toll–like Receptor 7 Agonist Imiquimod, in Patients on Dialysis. *Clin Infect Dis.* 2020 Jun 18;ciaa804. doi: 10.1093/cid/ciaa804. Online ahead of print.

CMV

Adler SP et al. Phase 1 Clinical Trial of a Conditionally Replication Defective Human Cytomegalovirus (CMV) Vaccine in CMV–Seronegative Subjects. *The Journal of Infectious Diseases* 2019;220:411–9. Read full article.

COVID

Fan S, et al. Ameliorated immunity elicited by intradermal inoculation in individuals vaccinated with inactivated SARS–CoV–2 vaccine. *Vaccine.* 2021 Nov 26;39(48):6980–6983. Read full article.

Levin Y. Peer review of: Immunogenicity and reactogenicity after booster dose with AZD1222 via intradermal route among adult who had received CoronaVac, Qeios, December 22, 2021. Read full article.

Kochba E. et al. Improved Insulin Pharmacokinetics Using a Novel Microneedle Device for Intradermal Delivery in Patients with Type 2 Diabetes. *Diabetes Technology & Therapeutics* 2016; 18(9). Read full article.

The background of the slide features a blue-toned image with a DNA double helix structure on the left and a network of interconnected nodes and lines on the right, suggesting a scientific or technological theme.

**THANKS FOR
WATCHING**

revival cell
wigomedi